



## SLMFN589 Series 589nm Single Longitudinal Mode Laser System up to 200mW

### Overview

The SLMFN589 series single longitudinal mode laser is ideal for applications that require a wavelength of 589nm and output power levels up to 200mW. The laser features narrow spectral linewidth, long coherent length, ultra-compact design, long lifetime, easy operating, and FDA-compliant system with driver. The laser is widely used in DNA sequencing, flow cytometry, cell sorting, optical instrument, spectrum analysis, interference, measurement, holography, physics experiment, and many other applications.



### Specifications

<b>Model Number</b>	<b>SLMFN589-XYP</b>		
Wavelength (nm)	589±1		
Operating mode	CW		
Output power (mW)	>30 (X=30), >50 (X=50), >100 (X=100)	>150 (X=150), >200 (X=200)	
Power stability (rms, over 4 hours)	<5% (P=A), <3% (P=E)	<5% (P=A), <3% (P=E), <2% (P=2)	
Transverse mode	TEM <sub>00</sub>		
Longitudinal mode	Single		
Spectral linewidth (nm)	<0.00001		
Coherent length (m)	>50		
Noise of amplitude (rms, 1Hz-20MHz)	<1%, typical <0.5%		
Beam diameter at aperture (1/e <sup>2</sup> , mm)	<2.5		
Beam divergence, full angle (mrad)	<1.2		
M <sup>2</sup> factor	<1.2		
Polarization ratio	>100:1, Horizontal ±5 degree		
Warm-up time (minutes)	<10		
Pointing stability after warm-up (mrad)	<0.05		
Beam height from base plate (mm)	27.4		
Operating temperature (°C)	15-35		
Dimensions of laser head (mm)	197(L)×70(W) ×50(H) mm <sup>3</sup>		
Weight of laser head (kg)	2.0 kg		
Power supply			
High Power Elite Power Supply (Y=H)	Input voltage	90-264VAC	
	Dimensions	275(L) ×145(W) ×104(H) mm <sup>3</sup>	
	Weight	2.3 kg	
	Notes	Fixed output power	
Expected lifetime (hours)	10,000		
Warranty period	10 months		
FDA Compliance	FDA CDRH Title 21 CFR 1040.10/11 Class IIIb		




#### Remarks:

- The laser head needs to be used on a heat sink with good heat dissipation.
- Specifications of the CW laser is based on the laser performance at full power output after the specified warmup period. The stability of output power may change when output power is adjusted at a different power level.






**Outline Dimensions (unit: mm)**

**Laser Head**



**Power Supply**



Note: The above specifications are subject to change without notice.